

EFHD1 Rabbit Polyclonal Antibody

Catalog	TDY069C	TDY069F	Tel: 010-82908854
Quantity	50 μ L	100 μ L	Free: 400-0620-621
			Web: www.tdybio.com

For research use only.

Applications	Species Cross-Reactivity	Molecular Weight	Isotype
WB	H	~27KD	IgG


Storage Buffer & Condition: PBS, pH 7.4, containing 0.02% **sodium azide** as Preservative and 50% Glycerol.
Store at **-20°C**. **Do not aliquot the antibody.**

Recommended dilutions: WB: 1:,000

Optimal dilutions should be determined by the end user.

Specificity: The EFHD1 Rabbit Polyclonal antibody detects endogenous EFHD1 proteins.

Background: EFHD1 is an EF-hand domain-containing protein that displays increased expression during neuronal differentiation.



Western blot analysis of 293T cell lysate with EFHD1 rabbit polyclonal antibody diluted at 1:2,000.