



## CD1 Mouse Monoclonal Antibody(9H6)

Catalog	TDY093C	TDY093F	Tel: 010-82908854
Quantity	50 $\mu$ L	100 $\mu$ L	Free: 400-0620-621
			Web: www.tdybio.com

*For research use only.*

Applications	Species	Cross-Reactivity	Molecular Weight	Isotype
IHC		H, R, M	N/A	IgG1

**Storage Buffer & Condition:** PBS, pH 7.4, containing 0.02% **sodium azide** as Preservative and 50% Glycerol.  
Store at **-20°C**. **Do not aliquot the antibody.**

**Recommended dilutions:** IHC: 1:200

**Optimal dilutions should be determined by the end user.**

**Specificity:** CD1 Mouse Monoclonal Antibody detects endogenous CD1 proteins.

**Background:** CD1 (cluster of differentiation 1) is a family of glycoproteins expressed on the surface of various human antigen-presenting cells. They are related to the class I MHC molecules, and are involved in the presentation of lipid antigens to T cells.