



Insulin Glargine Mouse Monoclonal Antibody(3F12)

Catalog TDY144C TDY144F

Tel: 010-82908854

Quantity 50 μ L 100 μ L

Free: 400-0620-621

Web: www.tdybio.com

For research use only.

Applications	Species	Cross-Reactivity	Molecular Weight	Isotype
ELISA		N/A	N/A	IgG1

Storage Buffer & Condition: PBS, pH 7.4, containing 0.02% **sodium azide** as Preservative and 50% Glycerol.

Store at **-20°C**. **Do not aliquot the antibody.**

Recommended dilutions: ELISA: 1:2,000-10,000

Optimal dilutions should be determined by the end user.

Specificity: Antibody can detect Insulin Glargine protein.

Alternative Names: N/A

Background: Insulin glargine, marketed by Sanofi-Aventis under the name Lantus, is a long-acting basal insulin analogue, given once daily to help control the blood sugar level of those with diabetes.