

## CBP-Tag Mouse Monoclonal Antibody( Mix-mA™)

Catalog	TDY154C	TDY154F	Tel: 010-82908854
Quantity	50μL	100μL	Free: 400-0620-621
			Web: www.tdybio.com

**For research use only.**

Applications	Species Cross-Reactivity	Molecular Weight	Isotype
WB, IP	N/A	N/A	IgG1

**Storage Buffer & Condition:** PBS, pH 7.4, containing 0.02% **sodium azide** as Preservative and 50% Glycerol.  
Store at **-20°C**. **Do not aliquot the antibody.**


**Recommended dilutions:** WB: 1:5,000 IP: 1:200

**Optimal dilutions should be determined by the end user.**

**Specificity:** Antibody can detects CBP recombinant protein.

**Alternative Names:** calmodulin binding peptide

**Background:** Epitope mapping within the calmodulin binding peptide (CBP) region of the TAP (tandem affinity purification) tag of human.



Western blot analysis of CBP recombinant protein with TDY154 diluted at

1) 1:5,000

2) 1:10,000