

ECFP Rabbit Polyclonal Antibody

Catalog TDY174C TDY174F

Tel: 010-82908854

Quantity 50 μ L 100 μ L

Free: 400-0620-621

Web: www.tdybio.com

For research use only.

Applications	Species Cross-Reactivity	Molecular Weight	Isotype
WB, IP	N/A	N/A	IgG

Storage Buffer & Condition: PBS, pH 7.4, containing 0.02% **sodium azide** as Preservative and 50% Glycerol.

Store at **-20°C**. **Do not aliquot the antibody.**

Recommended dilutions: WB: 1:2,000-5,000 IP: 1:200

Optimal dilutions should be determined by the end user.

Specificity: The ECFP Antibody can detects ECFP & CFP tag protein..

Alternative Names: Cyan fluorescent protein, Enhanced Cyan fluorescent protein

Background: Cyan Fluorescent Protein (CFP) is a genetic mutant of green fluorescent protein (GFP) originally derived from the jellyfish *Aequorea victoria*.



Western blot analysis of ECFP recombinant protein with TDY174 diluted at

1:5,000.