

TLR2 Rabbit Polyclonal Antibody(M25)

Catalog	TDY367C	TDY367F	Tel: 010-82908854
Quantity	50 μ L	100 μ L	Web: www.tdybio.com
			Entrez-Gene ID:7097 , Swiss-Prot Acc. O60603

For research use only.

Applications	Species Cross-Reactivity	Molecular Weight	Isotype
IHC	H, R, M	N/A	IgG

Storage Buffer & Condition: Antigen Affinity Purified IgG in PBS, pH 7.4, containing 0.02% **sodium azide** as Preservative and 50% Glycerol.

Store at **-20°C**. **Do not aliquot the antibody.**


Recommended dilutions: IHC: 1:200-500

Optimal dilutions should be determined by the end user.

Specificity: Antibody can detects endogenous TLR2 protein.

Alternative Names: CD282 antibody, TIL4 antibody, Toll like receptor antibody

Background: Toll-like receptor 2 also known as TLR2 is a protein that in humans is encoded by the *TLR2* gene. TLR2 has also been designated as CD282. TLR2 is one of the toll-like receptors and plays a role in the immune system.



Immunohistochemical analysis of paraffin-embedded Mouse Spleen Tissue using TLR2 (TDY367) Rabbit pAb diluted at 1:500.