

GAL4 Activation Domain Rabbit Polyclonal Antibody(AD)

Catalog TDY386C TDY386F

Tel: 010-82908854

Quantity 50 μ L 100 μ L

Web: www.tdybio.com

Entrez-Gene ID , Swiss-Prot Acc.

For research use only.

Applications	Species Cross-Reactivity	Molecular Weight	Isotype
WB	N/A	N/A	IgG

Storage Buffer & Condition: PBS, pH 7.4, containing 0.02% **sodium azide** as Preservative and 50% Glycerol.
Store at **-20°C**. **Do not aliquot the antibody.**


Recommended dilutions: WB: 1:10,000-50,000

Optimal dilutions should be determined by the end user.

Specificity: Antibody can detects recombinant GAL4 Activation Domain protein.

Alternative Names: N/A

Background: GAL4 is a transcriptional activator which controls the expression of genes involved in the catabolism of galactose.



Western blot analysis of Recombinant AD protein with GAL4 Activation

Domain Rabbit pAb TDY386 diluted at 1:20,000.