

Bcl-XI Rabbit Polyclonal Antibody(F267)

Catalog	TDY460C	TDY460F	Tel: 010-82908854
			Web: www.tdybio.com
Quantity	50μL	100μL	Entrez-Gene ID#598 , Swiss-Prot Acc.#Q07817

For research use only.

Applications	Species Cross-Reactivity	Molecular Weight	Isotype
IHC	H, R, M	30KD	IgG

Storage Buffer & Condition: Antigen Affinity Purified IgG in PBS, pH 7.4, containing 0.02% **sodium azide** as Preservative and 50% Glycerol.

Store at **-20°C**. **Do not aliquot the antibody.**


Recommended dilutions: IHC: 1:50-100

Optimal dilutions should be determined by the end user.

Specificity: Antibody can detects endogenous Bcl-XI protein.

Alternative Names: B2CL1 antibody, Bcl 2 like 1 antibody, Apoptosis regulator Bcl X antibody

Background: Bcl-xL prevents apoptosis through two different mechanisms: heterodimerization with an apoptotic protein inhibits its apoptotic effect and formation of mitochondrial outer membrane pores help maintain a normal membrane state under stressful conditions.



Immunohistochemical analysis of paraffin-embedded Human Lung Carcinoma using Bcl-XI(TDY460) Rabbit pAb diluted at 1:200.