

Cyclophilin B Mouse Monoclonal Antibody(2B10)

Catalog TDY487C TDY487F

Tel: 010-80117836

Web: www.tdybio.com

Quantity 50μL 100μL

Entrez-Gene ID#5479 , Swiss-Prot Acc.#P23284

For research use only.

Applications	Species Cross-Reactivity	Molecular Weight	Isotype
WB	H,R,M	21KD	IgG1

Storage Buffer & Condition: PBS, pH 7.4, containing 0.02% **sodium azide** as Preservative and 50% Glycerol.
Store at **-20°C. Do not aliquot the antibody.**


Recommended dilutions: WB: 1:1,000-2,000

Optimal dilutions should be determined by the end user.

Specificity: Antibody can detects endogenous Cyclophilin B protein.

Alternative Names: Cphn2, MGC2224, Rotamase, S cyclophilin, CYPB, CyP-S1, Rotamase B

Background: Cyclophilin B (CyPB) is an ER-localized chaperone protein belonging to the family of peptidyl-prolyl cis-trans isomerases.



Western blot analysis of 1) Jurkat, 2)293T, 3)Rat Liver Tissue, 4)3T3 with Cyclophilin B Mouse mAb TDY487 diluted at 1:2,000.