

## GRIK2(GluR6) Rabbit Polyclonal Antibody(A240)

Catalog TDY537C TDY537F

Tel: 010-80117836

Web: www.tdybio.com

Quantity 50μL 100μL

Entrez-Gene ID# 2898, Swiss-Prot Acc.#Q13002

**For research use only.**

Applications	Species Cross-Reactivity	Molecular Weight	Isotype
IHC	H	~103KD	IgG

**Storage Buffer & Condition:** Antigen Affinity Purified IgG in PBS, pH 7.4, containing 0.02% **sodium azide** as Preservative and 50% Glycerol.

Store at **-20°C**. **Do not aliquot the antibody.**


**Recommended dilutions:** IHC: 1:100-200

**Optimal dilutions should be determined by the end user.**

**Specificity:** Antibody can detects endogenous GRIK2(GluR6) protein.

**Alternative Names:** Excitatory amino acid receptor 4, GLR6, GLUK2, Glutamate receptor 6

**Background:** Ionotropic glutamate receptor. L-glutamate acts as an excitatory neurotransmitter at many synapses in the central nervous system.



Immunohistochemical analysis of paraffin-embedded Human Colon Tissue

using GRIK2(GluR6) (TDY537) Rabbit pAb diluted at 1:200.