



KChIP1 Rabbit Polyclonal Antibody(A249)

Catalog TDY544C TDY544F

Tel: 010-80117836

Web: www.tdybio.com

Quantity 50μL 100μL

Entrez-Gene ID# 30820, Swiss-Prot Acc.#Q9NZI2

For research use only.

Applications	Species Cross-Reactivity	Molecular Weight	Isotype
IHC	H,M	~30KD	IgG

Storage Buffer & Condition: Antigen Affinity Purified IgG in PBS, pH 7.4, containing 0.02% **sodium azide** as Preservative and 50% Glycerol.


Store at **-20°C. Do not aliquot the antibody.**

Recommended dilutions: IHC: 1:100-200

Optimal dilutions should be determined by the end user.

Specificity: Antibody can detects endogenous KChIP1 protein.

Alternative Names: KCNIP1,MGC9,MGC95,A type potassium channel modulatory protein 1



Immunohistochemical analysis of paraffin-embedded Mouse Brain Tissue
using KChIP1 (TDY544) Rabbit pAb diluted at 1:200.