



KChIP3 Rabbit Polyclonal Antibody(A250)

Catalog TDY545C TDY545F

Tel: 010-80117836

Web: www.tdybio.com

Quantity 50 μ L 100 μ L

Entrez-Gene ID# 30818, Swiss-Prot Acc.#Q9Y2W7

For research use only.

Applications	Species Cross-Reactivity	Molecular Weight	Isotype
IHC	H,R,M	~29KD	IgG

Storage Buffer & Condition: Antigen Affinity Purified IgG in PBS, pH 7.4, containing 0.02% **sodium azide** as Preservative and 50% Glycerol.


Store at **-20°C**. **Do not aliquot the antibody.**

Recommended dilutions: IHC: 1:100-200


Optimal dilutions should be determined by the end user.

Specificity: Antibody can detects endogenous KChIP3 protein.

Alternative Names: KCNIP3, A type potassium channel modulatory protein 3



Immunohistochemical analysis of paraffin-embedded Rat BrainTissue using KChIP3 (TDY545) Rabbit pAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Mouse BrainTissue using KChIP3 (TDY545) Rabbit pAb diluted at 1:200.