

Kir4.1 Rabbit Polyclonal Antibody(A252)

Catalog TDY546C TDY546F

Tel: 010-80117836

Web: www.tdybio.com

Quantity 50μL 100μL

Entrez-Gene ID# 3766, Swiss-Prot Acc.#P78508

For research use only.

Applications	Species Cross-Reactivity	Molecular Weight	Isotype
IHC	H,R,M	~42KD	IgG

Storage Buffer & Condition: Antigen Affinity Purified IgG in PBS, pH 7.4, containing 0.02% **sodium azide** as Preservative and 50% Glycerol.


Store at **-20°C. Do not aliquot the antibody.**

Recommended dilutions: IHC: 1:100-200


Optimal dilutions should be determined by the end user.

Specificity: Antibody can detects endogenous Kir4.1 protein.

Alternative Names: BIRK10,KCNJ10,potassium channel



Immunohistochemical analysis of paraffin-embedded Rat Kidney Tissue using Kir4.1 (TDY546) Rabbit pAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Mouse Brain Tissue using Kir4.1 (TDY546) Rabbit pAb diluted at 1:200.