



KV1.5 Rabbit Polyclonal Antibody(A257)

Catalog TDY549C TDY549F

Tel: 010-80117836

Web: www.tdybio.com

Quantity 50μL 100μL

Entrez-Gene ID# 3741, Swiss-Prot Acc.#P22460

For research use only.

Applications	Species Cross-Reactivity	Molecular Weight	Isotype
WB	H	~67KD	IgG

Storage Buffer & Condition: Antigen Affinity Purified IgG in PBS, pH 7.4, containing 0.02% **sodium azide** as Preservative and 50% Glycerol.


Store at **-20°C**. **Do not aliquot the antibody.**

Recommended dilutions: WB: 1:1,000-2,000

Optimal dilutions should be determined by the end user.

Specificity: Antibody can detects endogenous KV1.5 protein.

Alternative Names: KCNA5,PCN5,Potassium channel 1, voltage gated potassium channel protein kv1.5



Western blot analysis of 293T with KV1.5 Rabbit pAb TDY549 diluted at 1:2,000.