



KV10.1 Rabbit Polyclonal Antibody(A259)

Catalog TDY551C TDY551F

Tel: 010-80117836

Web: www.tdybio.com

Quantity 50 μ L 100 μ L

Entrez-Gene ID# 170850, Swiss-Prot Acc.#Q8TAE7

For research use only.

Applications	Species Cross-Reactivity	Molecular Weight	Isotype
IHC	H,R,M	~50KD	IgG

Storage Buffer & Condition: Antigen Affinity Purified IgG in PBS, pH 7.4, containing 0.02% **sodium azide** as Preservative and 50% Glycerol.


Store at **-20°C**. **Do not aliquot the antibody.**

Recommended dilutions: IHC: 1:100-200

Optimal dilutions should be determined by the end user.

Specificity: Antibody can detects endogenous KV10.1 protein.

Alternative Names: KCNG3,KV6.3,Voltage-gated potassium channel subunit kv10.1



Immunohistochemical analysis of paraffin-embedded Rat BrainTissue using Kv10.1 (TDY551) Rabbit pAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Mouse BrainTissue using Kv10.1 (TDY551) Rabbit pAb diluted at 1:200.