

## DUSP6 Rabbit Polyclonal Antibody

Catalog	TDY555C	TDY555F	Tel: 010-80117836
Quantity	50 $\mu$ L	100 $\mu$ L	Web: www.tdybio.com
			Entrez-Gene ID# 1848, Swiss-Prot Acc.#Q16828

**For research use only.**

Applications	Species Cross-Reactivity	Molecular Weight	Isotype
WB	H,R,M	~42,44KD	IgG

**Storage Buffer & Condition:** Antigen Affinity Purified IgG in PBS, pH 7.4, containing 0.02% **sodium azide** as Preservative and 50% Glycerol.


Store at **-20°C. Do not aliquot the antibody.**

**Recommended dilutions: WB: 1:1,000-5,000**

**Optimal dilutions should be determined by the end user.**

**Specificity:** Antibody can detects endogenous DUSP6 protein.

**Alternative Names:** Dual specificity phosphatase 6, HH19, MKP3



Western blot analysis of 1) HepG2, 2)3T3 , 3) Rat Heart Tissue with DUSP6

Rabbit pAb TDY555 diluted at 1:3,000.