

## PPAR Delta Mouse Monoclonal Antibody(1D7)

Catalog TDY571C TDY571F

Tel: 010-80117836

Web: [www.tdybio.com](http://www.tdybio.com)

Quantity 50μL 100μL

Entrez-Gene ID#5469 , Swiss-Prot Acc.#Q03181

**For research use only.**

| Applications | Species Cross-Reactivity | Molecular Weight | Isotype |
|--------------|--------------------------|------------------|---------|
| IHC          | H,R,M                    | ~50KD            | IgG1    |

**Storage Buffer & Condition:** PBS, pH 7.4, containing 0.02% **sodium azide** as Preservative and 50% Glycerol.

Store at **-20°C**. **Do not aliquot the antibody.**


**Recommended dilutions:** IHC: 1:100-200

**Optimal dilutions should be determined by the end user.**


**Specificity:** The PPAR Delta antibody can detects endogenous PPAR Delta protein.

**Alternative Names:** FAAR, NR1C2, NUC1, Peroxisome proliferative activated receptor delta

**Background:**



Immunohistochemical analysis of paraffin-embedded Human Brain Tissue using PPAR Delta (TDY571) Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Mouse Brain Tissue using PPAR Delta (TDY571) Mouse mAb diluted at 1:200.