

Cystatin C Mouse Monoclonal Antibody(7F11)

Catalog TDY594C TDY594F

Tel: 010-80117836

Web: www.tdybio.com

Quantity 50μL 100μL

Entrez-Gene ID#1471 , Swiss-Prot Acc.#P01034

For research use only.

Applications	Species Cross-Reactivity	Molecular Weight	Isotype
WB,IHC,ELISA	H	~14KD	IgG1

Storage Buffer & Condition: PBS, pH 7.4, containing 0.02% **sodium azide** as Preservative and 50% Glycerol.


Store at **-20°C**. **Do not aliquot the antibody.**

Recommended dilutions: WB: 1:1,000-2,000 IHC: 1:100-200 ELISA: 1:10,000-20,000


Optimal dilutions should be determined by the end user.

Specificity: Antibody can detects endogenous Cystatin C protein.

Alternative Names: AD8, Cystatin3, Gamma trace, HCCAA



Immunohistochemical analysis of paraffin-embedded Human Kidney Tissue using Cystatin C (TDY594) Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Brain Tissue using Cystatin C (TDY594) Mouse mAb diluted at 1:200.