



Cystatin C Mouse Monoclonal Antibody(3C6)

Catalog	TDY596C	TDY596F	Tel: 010-80117836
			Web: www.tdybio.com
Quantity	50μL	100μL	Entrez-Gene ID#1471 , Swiss-Prot Acc.#P01034

For research use only.

Applications	Species Cross-Reactivity	Molecular Weight	Isotype
WB,ELISA	H	~14KD	IgG1

Storage Buffer & Condition: PBS, pH 7.4, containing 0.02% **sodium azide** as Preservative and 50% Glycerol.
Store at **-20°C**. **Do not aliquot the antibody.**

Recommended dilutions: WB: 1:1,000-2,000 ELISA: 1:10,000-20,000

Optimal dilutions should be determined by the end user.

Specificity: Antibody can detects endogenous Cystatin C protein.

Alternative Names: AD8, Cystatin3, Gamma trace, HCCAA