

Muscle Actin Mouse Monoclonal Antibody(1E9)

Catalog TDY605C TDY605F

Tel: 010-80117836

Web: www.tdybio.com

Quantity 50μL 100μL

Entrez-Gene ID#PAN , Swiss-Prot Acc.#PAN

For research use only.

Applications	Species Cross-Reactivity	Molecular Weight	Isotype
IHC	H,R,M	~42KD	IgG1

Storage Buffer & Condition: PBS, pH 7.4, containing 0.02% **sodium azide** as Preservative and 50% Glycerol.
Store at **-20°C**. **Do not aliquot the antibody.**


Recommended dilutions: IHC: 1:100-200

Optimal dilutions should be determined by the end user.


Specificity: The Muscle Actin antibody can detects endogenous Muscle Actin protein.

Alternative Names:


Background:



Immunohistochemical analysis of paraffin-embedded Human Colon Carcinoma Tissue using Muscle Actin (TDY605) Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Colon Carcinoma Tissue using Muscle Actin (TDY605) Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Mouse Skeletal Muscle Tissue using Muscle Actin (TDY605) Mouse mAb diluted at 1:200.

Applications: WB-Western blot IHC-Immunochemistry IF-Immunofluorescence IP-Immunoprecipitation ChIP-Chormatin Immunoprecipitation
Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish