



Luciferase Mouse Monoclonal Antibody(6B8)

Catalog	TDY614C	TDY614F	Tel: 010-80117836
			Web: www.tdybio.com
Quantity	50μL	100μL	Entrez-Gene ID# , Swiss-Prot Acc. #P08659

For research use only.

Applications	Species Cross-Reactivity	Molecular Weight	Isotype
WB	Firefly	~60KD	IgG1

Storage Buffer & Condition: PBS, pH 7.4, containing 0.02% **sodium azide** as Preservative and 50% Glycerol.
Store at **-20°C**. **Do not aliquot the antibody.**

Recommended dilutions: WB: 1:1,000-2,000

Optimal dilutions should be determined by the end user.

Specificity: The Luciferase antibody can detect endogenous and recombinant Luciferase protein.

Alternative Names:

Background: