

Epsilon Tubulin Mouse Monoclonal Antibody(2B3)

Catalog TDY622C TDY622F

Tel: 010-80117836

Web: www.tdybio.com

Quantity 50μL 100μL

Entrez-Gene ID#51175 , Swiss-Prot Acc.#Q9UJT0

For research use only.

Applications	Species Cross-Reactivity	Molecular Weight	Isotype
WB	H,R,M	~55KD	IgG1

Storage Buffer & Condition: PBS, pH 7.4, containing 0.02% **sodium azide** as Preservative and 50% Glycerol.

Store at **-20°C. Do not aliquot the antibody.**


Recommended dilutions: WB: 1:1,000-2,000

Optimal dilutions should be determined by the end user.

Specificity: The epsilon tubulin antibody can detects endogenous epsilon tubulin protein.

Alternative Names: FLJ22589,TUBE1,tubulin epsilon 1, ε tubulin

Background:



Western blot analysis of 1) MCF7 Cell Lysate, 2) Rat Brain Tissue Lysate, 3)

Mouse Brain Tissue Lysate using Epsilon Tubulin (TDY622) Mouse mAb

diluted at 1:2000.