



Luciferase Mouse Monoclonal Antibody(2G11)

Catalog TDY624C TDY624F

Tel: 010-80117836

Web: www.tdybio.com

Quantity 50 μ L 100 μ L

Entrez-Gene ID# , Swiss-Prot Acc.#P08659

For research use only.

Applications	Species Cross-Reactivity	Molecular Weight	Isotype
WB	Firefly	~60KD	IgG1

Storage Buffer & Condition: PBS, pH 7.4, containing 0.02% **sodium azide** as Preservative and 50% Glycerol.

Store at **-20°C**. **Do not aliquot the antibody.**


Recommended dilutions: WB: 1:1,000-2,000

Optimal dilutions should be determined by the end user.

Specificity: The Luciferase antibody can detects endogenous and recombinant Luciferase protein.

Alternative Names:

Background:



Western blot analysis of Firefly Luciferase Transfected HeLa Cell Lysate
using Luciferase (TDY624) Mouse mAb diluted at 1:2000.