

Epsilon Tubulin Mouse Monoclonal Antibody(3C1)

Catalog TDY630C TDY630F

Tel: 010-80117836

Web: www.tdybio.com

Quantity 50μL 100μL

Entrez-Gene ID#51175 , Swiss-Prot Acc.#Q9UJT0

For research use only.

Applications	Species Cross-Reactivity	Molecular Weight	Isotype
IHC	H,R,M	~55KD	IgG1

Storage Buffer & Condition: PBS, pH 7.4, containing 0.02% **sodium azide** as Preservative and 50% Glycerol.
Store at **-20°C**. **Do not aliquot the antibody.**


Recommended dilutions: IHC: 1:100-200

Optimal dilutions should be determined by the end user.


Specificity: The epsilon tubulin antibody can detects endogenous epsilon tubulin protein.

Alternative Names: FLJ22589,TUBE1,tubulin epsilon 1, ε tubulin

Background:



Immunohistochemical analysis of paraffin-embedded Human Colon Carcinoma Tissue using Epsilon Tubulin (TDY630) Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma Tissue using Epsilon Tubulin (TDY630) Mouse mAb diluted at 1:200.