

JAK1 Rabbit Polyclonal Antibody(F172)

Catalog TDY687C TDY687F

Tel: 010-80117836

Web: www.tdybio.com

Quantity 50μL 100μL

Entrez-Gene ID#3716 , Swiss-Prot Acc.#P23458

For research use only.

Applications	Species Cross-Reactivity	Molecular Weight	Isotype
IHC	H,R,M	~130KD	IgG

Storage Buffer & Condition: Antigen Affinity Purified IgG in PBS, pH 7.4, containing 0.02% **sodium azide** as Preservative and 50% Glycerol.

Store at **-20°C. Do not aliquot the antibody.**



Recommended dilutions: IHC: 1:100-200

Optimal dilutions should be determined by the end user.

Specificity: The Antibody can detects endogenous JAK1 proteins.

Alternative Names: Janus kinase 1,JTK3,Tyrosine protein kinase JAK1,JAK1A,JAK1B

Background: Members of the Janus family of tyrosine kinases (Jak1, Jak2, Jak3 and Tyk2) are activated by ligands binding to a number of associated cytokine receptors.



Immunohistochemical analysis of paraffin-embedded Human Lung Carcinoma Tissue using JAK1 (TDY687) Rabbit pAb diluted at 1:200

Immunohistochemical analysis of paraffin-embedded Human Lung Carcinoma Tissue using JAK1 (TDY687) Rabbit pAb diluted at 1:200