

mScarlet Mouse Monoclonal Antibody(Mix-mA)

Catalog	TDY1079C	TDY1079F	Tel: 010-80117836
			Web: www.tdybio.com
Quantity	50μL	100μL	Entrez-Gene ID# , Swiss-Prot Acc.#

For research use only.

Applications	Species Cross-Reactivity	Molecular Weight	Isotype
WB	Recombinant		IgG1

Storage Buffer & Condition: PBS, pH 7.4, containing 0.02% **sodium azide** as Preservative and 50% Glycerol.
Store at **-20°C. Do not aliquot the antibody.**


Recommended dilutions:WB: 1:2,000-5,000

Optimal dilutions should be determined by the end user.

Specificity:Antibody can detects recombinant mScarlet proteins.

Alternative Names:

Background: mScarlet: a bright monomeric red fluorescent protein for cellular imaging.



Western blot analysis of mScarlet recombinant transfected 293T Cell Lysate
using mScarlet (TDY1079) Mouse mAb diluted at 1:5,000.