

mScarlet Rabbit Polyclonal Antibody

Catalog	TDY1087C	TDY1087F	Tel: 010-80117836
			Web: www.tdybio.com
Quantity	50μL	100μL	Entrez-Gene ID# , Swiss-Prot Acc.#

For research use only.

Applications	Species Cross-Reactivity	Molecular Weight	Isotype
WB	Recombinant		IgG

Storage Buffer & Condition: PBS, pH 7.4, containing 0.02% **sodium azide** as Preservative and 50% Glycerol.
Store at **-20°C. Do not aliquot the antibody.**


Recommended dilutions: WB: 1:1,000-2,000

Optimal dilutions should be determined by the end user.

Specificity: Antibody can detect recombinant mScarlet proteins.

Alternative Names:

Background: mScarlet: a bright monomeric red fluorescent protein for cellular imaging.



Western blot analysis of mScarlet recombinant transfected 293T Cell Lysate
using mScarlet (TDY1087) Rabbit Polyclonal Ab diluted at 1:2,000.