

JAK2 Mouse Monoclonal Antibody(6D6)

Catalog	TDY1062C	TDY1062F	Tel: 010-80117836
			Web: www.tdybio.com
Quantity	50μL	100μL	Entrez-Gene ID:3717 , Swiss-Prot Acc.O60674

For research use only.

Applications	Species Cross-Reactivity	Molecular Weight	Isotype
IHC	H,R,M	~120KD	IgG1

Storage Buffer & Condition: PBS, pH 7.4, containing 0.02% **sodium azide** as Preservative and 50% Glycerol.
Store at **-20°C**. **Do not aliquot the antibody.**


Recommended dilutions: IHC: **1:100-200**

Optimal dilutions should be determined by the end user.

Specificity: The Antibody can detects endogenous JAK2 proteins.

Alternative Names: Janus kinase 2, JTK10, Tyrosine protein kinase Jak2, THCYT3

Background: Members of the Janus family of tyrosine kinases (Jak1, Jak2, Jak3 and Tyk2) are activated by ligands binding to a number of associated cytokine receptors.



Immunohistochemical analysis of paraffin-embedded Human Lung Carcinoma Tissue using JAK2 (TDY1062) Mouse Monoclonal antibody diluted at 1:200.