

Chicken Ovalbumin Rabbit Polyclonal Antibody

Catalog	TDY858C	TDY858F	Tel: 010-80117836
			Web: www.tdybio.com
Quantity	50μL	100μL	Entrez-Gene ID# 396058, Swiss-Prot Acc.#P01012

For research use only.

Applications	Species Cross-Reactivity	Molecular Weight	Isotype
WB, E	Chicken	43KD	IgG

Storage Buffer & Condition: Antigen Affinity Purified IgG in PBS, pH 7.4, containing 0.02% **sodium azide** as Preservative and 50% Glycerol.

Store at **-20°C. Do not aliquot the antibody.**

Recommended dilutions: WB: 1:10,000-50,000 ELISA 1:10,000-50,000

Optimal dilutions should be determined by the end user.

Specificity: Antibody can detect endogenous Chicken Ovalbumin protein.

Alternative Names: OVA, OVAL, SERPINB 14, Egg albumin

Background: Ovalbumin is a protein of unknown function found in large quantities in avian egg-white.



Western blot analysis of Chicken Ovalbumin protein using Ovalbumin

(TDY858) Rabbit pAb diluted at 1:20,000.