

Branchiostoma blecheri PAX9 Rabbit Polyclonal Antibody

Catalog TDY860C TDY860F

Tel: 010-80117836

Web: www.tdybio.com

Quantity 50μL 100μL

Entrez-Gene ID# , Swiss-Prot Acc.#

For research use only.

Applications	Species Cross-Reactivity	Molecular Weight	Isotype
WB, IF	Branchiostoma blecheri (H, M, R)	38KD	IgG

Storage Buffer & Condition: Antigen Affinity Purified IgG in PBS, pH 7.4, containing 0.02% **sodium azide** as Preservative and 50% Glycerol.

Store at **-20°C**. **Do not aliquot the antibody.**


Recommended dilutions: WB: 1:500-1,000 IF 1:100-200

Optimal dilutions should be determined by the end user.

Specificity: Antibody can detect endogenous PAX9 protein.

Alternative Names: Paired box gene 9, STHAG3, Paired box 9

Background: Transcription factor required for normal development of thymus, parathyroid glands, ultimobranchial bodies, teeth, skeletal elements of skull and larynx as well as distal limbs.



Western blot analysis of 1) PC-3 Cell, 2) HepG2 Cell, 3) Jurkat Cell Lysate using PAX9 (TDY860) Rabbit pAb diluted at 1:1,000.