

SV40 T-antigen Rabbit Polyclonal Antibody

Catalog TDY872C TDY872F

Tel: 010-80117836

Web: www.tdybio.com

Quantity 50μL 100μL

Entrez-Gene ID# , Swiss-Prot Acc.#

For research use only.

Applications	Species Cross-Reactivity	Molecular Weight	Isotype
WB		KD	IgG

Storage Buffer & Condition: PBS, pH 7.4, containing 0.02% **sodium azide** as Preservative and 50% Glycerol.

Store at **-20°C**. **Do not aliquot the antibody.**

Recommended dilutions: WB: 1:10,000-20,000

Optimal dilutions should be determined by the end user.

Specificity: Antibody can detect Recombinant SV40 T-antigen proteins.

Alternative Names:

Background:



Western blot analysis of Recombinant protein using SV40 T-antigen

(TDY872) Rabbit pAb diluted at 1:10,000.