

## Luciferase Rabbit Polyclonal Antibody(9D6)

|          |         |         |                                          |
|----------|---------|---------|------------------------------------------|
| Catalog  | TDY873C | TDY873F | Tel: 010-80117836                        |
|          |         |         | Web: www.tdybio.com                      |
| Quantity | 50μL    | 100μL   | Entrez-Gene ID# , Swiss-Prot Acc.#P08659 |

**For research use only.**

| Applications | Species Cross-Reactivity | Molecular Weight | Isotype |
|--------------|--------------------------|------------------|---------|
| WB           | Firefly                  | ~60KD            | IgG1    |

**Storage Buffer & Condition:** PBS, pH 7.4, containing 0.02% **sodium azide** as Preservative and 50% Glycerol.  
Store at **-20°C**. **Do not aliquot the antibody.**

**Recommended dilutions:** WB: 1:1,000-2,000

**Optimal dilutions should be determined by the end user.**

**Specificity:** The Luciferase antibody can detects endogenous and recombinant Luciferase protein.

**Alternative Names:**

**Background:**



Western blot analysis of Recombinant FireFly protein using Luciferase  
(TDY873) Rabbit pAb diluted at 1:2,000.