



SARS-COV-2 Nucleocapsid Protein Rabbit Polyclonal Antibody

Catalog TDY1107C TDY1107F

Quantity 50μL 100μL

Tel: 010-80117836

Web: www.tdybio.com

Entrez-Gene ID:N/A , Swiss-Prot Acc.

For research use only.

Applications	Species Cross-Reactivity	Molecular Weight	Isotype
WB	Species independent		IgG

Storage Buffer & Condition: PBS, pH 7.4, containing 0.02% **sodium azide** as Preservative and 50% Glycerol.
Store at **-20°C**. **Do not aliquot the antibody.**

Recommended dilutions: WB: 1:5,000-10,000

Optimal dilutions should be determined by the end user.

Specificity: The Antibody can detect SARS-COV-2 Nucleocapsid proteins.

Alternative Names: N protein, 2019-nCoV nucleocapsid protein

Background: Nucleocapsid protein is a most abundant protein of coronavirus on the helical nucleocapsid of coronaviruses.