

Syndecan-4 Rabbit Polyclonal Antibody

Catalog TDY1110C TDY1110F

Quantity 50μL 100μL

Tel: 010-80117836

Web: www.tdybio.com

Entrez-Gene ID:6385 , Swiss-Prot Acc.P31431

For research use only.

Applications	Species Cross-Reactivity	Molecular Weight	Isotype
WB	H, M, R	26KD	IgG

Storage Buffer & Condition: PBS, pH 7.4, containing 0.02% **sodium azide** as Preservative and 50% Glycerol.
Store at **-20°C. Do not aliquot the antibody.**


Recommended dilutions:WB: 1:1,000-2,000

Optimal dilutions should be determined by the end user.

Specificity:The Antibody can detect endogenous Syndecan-4 proteins.

Alternative Names: SDC4; YND4; Amphiglycan; Ryudocan core protein

Background: The protein encoded by this gene is a transmembrane (type I) heparan sulfate proteoglycan that functions as a receptor in intracellular signaling. The encoded protein is found as a homodimer and is a member of the syndecan proteoglycan family.



Western blot analysis of Mouse Liver Tissue Lysate using Syndecan-4 (TDY1110) Rabbit Polyclonal Antibody diluted at 1:2,000.